

ALABAMA DEPARTMENT OF FORENSIC SCIENCES TOXICOLOGY DISCIPLINE

Biological Specimens Kits Specifications

Antemortem specimens (DUI) collection kits


- 1 ea Kit box, white corrugated cardboard, approx 5" x 6.5" x 2.5", with affixed integrity seal; to include on address face as stamped or affixed label: "ALABAMA DEPARTMENT OF FORENSIC SCIENCES DUI / ANTEMORTEM SPECIMENS COLLECTION KIT"; kit box must comply with common carrier and US Postal regulations pertaining to shipment of biological specimens; kit box to include on address face as stamped or affixed label: "EXEMPT HUMAN SPECIMENS"
- 1 ea Kit box to include as stamped or affixed label: "THIS KIT PURCHASED THROUGH THE ALABAMA CHEMICAL TESTING, TRAINING AND EQUIPMENT TRUST FUND", vendor's kit part/reorder and lot numbers and expiration date; kit box label to state that the kit expiration applies to blood collection tubes and that expired tubes may be replaced with current hospital or clinical stock
- 1 ea FDA insert
- 1 ea Instruction sheet, 1-page 8.5" x 11", verbiage to be provided by purchaser
- 1 ea Cleansing pad, alcohol free
- 2 ea Blood collection tube, glass, vacuum, gray-stopper, 10 mL, sterile, containing 100 mg of sodium fluoride and 20 mg of potassium oxalate, with affixed label to provide space for hand-written identification; affixed label also to identify contents of empty tube (typical manufacturer's label is sufficient)
- 1 ea Specimen bottle, non-opaque plastic, screw-cap, 100 mL, with affixed integrity seal, with affixed label to provide space for hand-written identification
- 1 ea Holder/protector insert(s), to separate and firmly secure all specified specimen containers within the kit box
- 1 ea Hypodermic needle, 20 gauge x 1.5", sterile
- 1 ea Plastic bag, sealing (Ziplok® or otherwise), 8" x 8" or sufficient to enclose holder/protector with filled specimen containers, to include liquid-absorbing sheet; bag to include international biohazard symbol as stamped, printed or affixed label
- 3 ea Specimen security seal, 1" x 7" approx, tamper-evident, white to provide for hand-written labeling; per attachment (example)
- 1 ea Kit shipping seal, white, tamper-evident, for re-sealing kit after specimen collection

Vendor to provide a representative example of the proposed kit

ALABAMA DEPARTMENT OF FORENSIC SCIENCES TOXICOLOGY DISCIPLINE

Biological Specimens Kits Specifications

Antemortem specimens (DUI) collection kits – Container seals

<hr/> Subject's Name <hr/>		SPECIMEN SECURITY SEAL	<hr/> Investigating Officer's Initials <hr/>
<hr/> Collector's Initials <hr/>			<hr/> Agency/Police Dept. <hr/>
<hr/> Date <hr/>	<hr/> Time <hr/>	BIOHAZARD (PLACE OVER BOTTLE CAP/BLOOD TUBE STOPPER)	<hr/> Identification No. #165A <hr/>
<hr/> Subject's Name <hr/>		SPECIMEN SECURITY SEAL	<hr/> Investigating Officer's Initials <hr/>
<hr/> Collector's Initials <hr/>			<hr/> Agency/Police Dept. <hr/>
<hr/> Date <hr/>	<hr/> Time <hr/>	BIOHAZARD (PLACE OVER BOTTLE CAP/BLOOD TUBE STOPPER)	<hr/> Identification No. #165A <hr/>
<hr/> Subject's Name <hr/>		SPECIMEN SECURITY SEAL	<hr/> Investigating Officer's Initials <hr/>
<hr/> Collector's Initials <hr/>			<hr/> Agency/Police Dept. <hr/>
<hr/> Date <hr/>	<hr/> Time <hr/>	BIOHAZARD (PLACE OVER BOTTLE CAP/BLOOD TUBE STOPPER)	<hr/> Identification No. #165A <hr/>